New Cell-Based High-Throughput Screening Platform for RNA-Mediated Neurodegenerative Diseases

(In collaboration with Dr. Samie R. Jaffrey, Weill Medical College of Cornell University)

Introduction

Most neurodegenerative disorders are thought to be caused by misfolded proteins that form toxic aggregates and disrupt cellular homeostasis in neurons. However, recent advances have identified abnormal expansion of trinucleotide repeat sequences in disease genes as a new RNA-mediated mechanism that plays a causative role in neurodegeneration. Specifically, toxic expanded RNA sequences form large nuclear "sinks" to sequester RNA-binding proteins and splicing factors that prevent normal cellular functions and induce neuronal cell death. This RNA-mediated mechanism is now associated with at least 10 neurological disorders that include several forms of inherited ataxias, Huntington disease, and myotonic dystrophies. Currently, there is no drug approved to slow the progression or reverse the pathology for many if not most of the RNA-mediated diseases. The majority of the medically available treatments only serve to manage symptoms and disability. Therefore, there is a strong need to discover first-in-class drug treatments capable of preventing and/or reversing these disease progression.

Our Solution

A new screening platform to accelerate the drug discovery of RNA-mediated neurodegenerative disease:


- A cell-based high-throughput screening system that will enable rapid identification of drugs that can specifically disrupt toxic RNA formations.
- A drug discovery platform that is both highly selective for the toxic RNA of interest and can be broadly applied across most RNA-mediated diseases.
- Compatible with standard fluorescence-based microtiter plate readers found in most HTS facilities

We believe that our technology offers an unique access to a new class of therapeutic targets and addresses a significant unmet medical market.

Karen Y. Wu

Lucerna, Inc., New York, NY

Results






Figure 3. Live-cell imaging of RNA foci development in human fibroblasts. (A) Time-lapsed images of a COS7 cell transiently transfected with a CGG_{mut}-Spinach2 vector. Time 0 indicates the first frame that displayed fluorescence above background. White arrowheads mark small foci formed de novo; red and blue asterisks mark emerging foci. (B) Images of a cell containing CGG_{mut}-Spinach2 aggregates after treatment with 1 µg/ml actinomycin D, a transcription inhibitor. Spinach2 signal was stable and remained unchanged for up to 8 h.

ALS/FTD

SCA3

SCA8

SCA10

SCA12

HDL2


GGGGCC

Diseases with RNA toxicity as a component

CAG

CUG

CUG



Karen Y. Wu, Ph.D. karen.wu@lucernatechnologies.com

Lucerna, Inc. 3960 Broadway Suite 330E New York, NY 10032

Table 1. A list of neurological diseases where expanded-repeat RNA toxicity plays an important role in disease pathology.

7-28

500-4.500

66-78

Muscle weakness, cogniti

Ataxia and parkinsonism

Ataxia and slurred speech

Ataxia and seizures

Chorea, cognitive deficits

Ataxia, tremor, and dementia

Conclusions

Demonstration of a new approach to identity drugs that target RNA-mediated neurodegenerative diseases:

We have developed a novel technology that enables real-time imaging RNA movements in living cells.

We demonstrate that we can monitor the formation and dissolution of toxic ribonuclear foci that contain expanded CGG transcripts and RNA-binding proteins.

• Our platform can report the effects of small molecule inhibitors on toxic RNA foci aggregation

• Our platform can be easily adapted for other RNAmediated neurodegenerative diseases including ones that are listed in Table 1.

Publications

Paige JS, Wu KY, Jaffrey SR. (2011) RNA mimic of green fluorescent protein. Science.

Strack RL, Disney MD, Jaffrey SR. (2013) A superfolding Spinach2 reveals the dynamic nature of trinucleotide repeat RNA. Nature Methods.

Strack RL, Jaffrey SR. (2014) Using RNA mimics of GFP to image RNA dynamics in mammalian cells. Advanced Fluorescence *Microscopy*, 1st ed. Academic Press, Elsevier, Inc.

Contact Information

Telephone: (855) 582-3762 www.lucernatechnologies.com